

Proficiency Tests 2023: Registration Form

Akkreditierung / Accreditation

akkreditiert nach / accredited according to: ISO/IEC 17.043-2010



Die Akkreditierung gilt nur für mit Hinweis gekennzeichnete Eignungsprüfungs-Programme des in der Urkundenanlage genannten Umfangs.

The accreditation is valid only for proficiency test programs marked with a note according to the scope of the annex to the certificate of accreditation

Please, mark the respective proficiency tests in the related column for registration! And send the form by email to pt@dla-lvu.de or register directly on our website form at www.dla-lvu.de!

Get 5 % discount per ≥ 10 registrations until 15th of January 2023

1. Allergens - Elicitors of Allergies and Intolerance Reactions - Accredited Programs -

Register	Art.-No. (Lot-No.)	Proficiency Test (PT) 2023	Shipping Samples	Price
Allergen-PTs (qualitative / quantitative with z-scores and recovery rates) - Accredited Programs - <i>ptAL01 - 09: Two Samples (25g each) with Allergens in a Food Matrix, as pastry, chocolate ice cream or boiled sausage are provided for qualitative and quantitative determination. One of these samples contains the incurred allergen(s). The third "Spiking Level Sample" (15g) contains the same allergens in simple matrix without further processing with similar amounts. The evaluation of results will be qualitative and quantitative. For information DLA calculates a recovery rate and recovery z-scores for the "positive"-sample and the "Spiking Level Sample". Any suitable method can be applied for analysis.</i>				
<input type="checkbox"/>	ptAL01 (Lot 23CrMiSo)	Allergens I: Crustaceae, Milk (Casein) and Soy in Instant Soup Powder with "Spiking Level Sample" (3 different Samples)	January	315 €
<input type="checkbox"/>	ptAL02 (Lot 23LuSeGl)	Allergens II: Lupine, Black Sesame and Wheat (Gluten) in "gluten-free" Rice Cracker (baked) with "Spiking Level Sample" (3 different Samples)	February	315 €
<input type="checkbox"/>	ptAL03 (Lot 23MiGl)	Allergens III: β-Lactoglobulin, Casein and Gluten in Infant Food with "Spiking Level Sample" (3 different Samples)	April	315 €
<input type="checkbox"/>	ptAL04 (Lot 23CeMuSe)	Allergens IV: Celery, Mustard and Sesame in Boiled Sausage with "Spiking Level Sample" (3 different Samples)	March	315 €
<input type="checkbox"/>	ptAL05 (Lot 23CwHaEg)	Allergens V: Cashew, Hazelnut and Egg in Cocoa Biscuit (baked) with "Spiking Level Sample" (3 different Samples)	September	315 €
<input type="checkbox"/>	ptAL06 (Lot 23PeMaWa)	Allergens VI: Peanut, Macadamia and Walnut in Ice cream (Chocolate) with "Spiking Level Sample" (3 different Samples)	October	315 €
<input type="checkbox"/>	ptAL07 (Lot 23PiFiCe)	Allergens VII: Pistachio, Fish and Celery in Instant Soup Powder with "Spiking Level Sample" (3 different Samples)	November	315 €
<input type="checkbox"/>	ptAL08 (Lot 23AlBrPeBu)	Allergens VIII: Almond, Brazil Nut, Pecan and Buckwheat in Cereal Muesli with "Spiking Level Sample" (3 different Samples)	December	315 €
<input type="checkbox"/>	ptAL09 (Lot 23PeAlMo)	Allergens IX: Peanut, Almond and Mollusks in Barbecue Spice Mix (with onion, garlic and paprika powder) with "Spiking Level Sample" (3 different Samples)	June	315 €
<input type="checkbox"/>	ptAL10 (Lot 23GlSo)	Allergens X: Gluten and Soya in Oat Drink and Rice Drink (4 different Samples)	November	315 €
<input type="checkbox"/>	ptAL11 (Lot 23GlMuSo)	Allergens XI: Gluten, Mustard and Soya in Oat Flour and Rice Flour (4 different Samples)	July	315 €

Subscribe	Art.-No. (Lot-No.)	Proficiency Test (PT) 2023	Shipping Samples	Price
Allergen-Screening (qualitative) <i>ptALS1-3: 4 Samples (20g each) with Allergens in a simple Carrier Matrix are provided. The evaluation of results is strictly qualitative with indication of the obtained agreements with the consensus values of the participants and the spiking of samples. For analysis ELISA (+ Lateral Flow), PCR and LC/MS can be applied. Quantitative spiking levels are given in the reports.</i>				
<input type="checkbox"/> ELISA <input type="checkbox"/> PCR <input type="checkbox"/> LC/MS	ptALS1 (Lot 23)	Allergen-Screening I - 4 Samples qualitative: Cashew, Hazelnut, Coconut, Macadamia, Almond, Brazil Nuts, Pecan, Pistachio, Walnut	March	380 € ¹ 460 € ²
<input type="checkbox"/> ELISA <input type="checkbox"/> PCR <input type="checkbox"/> LC/MS	ptALS2 (Lot 23)	Allergen-Screening II - 4 Samples qualitative: Crustaceae, Egg, Fish, Milk, Molluscs, Mustard (yellow/white, brown and black), Soybean	June	380 € ¹ 460 € ²
<input type="checkbox"/> ELISA <input type="checkbox"/> PCR <input type="checkbox"/> LC/MS	ptALS3 (Lot 23)	Allergen-Screening III - 4 Samples qualitative: Cereals containing Gluten (Barley or Oat, Wheat and Rye), Peanut, Lupine, Buckwheat, Celery (Leaves / Stem, Root and Seed), Sesame (white and black)	October	380 € ¹ 460 € ²
<i>prices 1) participation with one method, 2) participation with 2 or 3 methods</i>				
Allergen-Response of processed allergenic Ingredients <i>ptALR1-2: 5 Samples (20g each) of a differently processed Allergen are provided for qualitative and quantitative determination, a simple Carrier Matrix as well as a „Blank Sample“ (Carrier Matrix). The individual processed ingredients of the samples are mixed gravimetrically by DLA each from 4 industrially manufactured products. The analytical methods are optional.</i>				
<input type="checkbox"/>	ptALR1 (Lot 23Cr)	Response PT Crustaceae : Processed samples of Crustaceae (raw), Crustaceae (cooked), Shrimp (fried), Shrimp Spread and Shrimp paste (fermented) in potato powder matrix (levels: 25 - 150 mg/kg)	April	400 €
<input type="checkbox"/>	ptALR2 (Lot 23Fi)	Response PT Fish : Processed samples Fish (raw), Fish (cooked), Fish Burger (fried), Surimi and Fish Sauce (fermented) in potato powder matrix (levels: 25 - 150 mg/kg)	November	400 €
<i>More Allergen-Response Single-PT-Sets are available for method validation and verification! See also our SPT programme</i>				
Action Level Matrix- ALM Verification <i>ptALM1-2: A kind of Calibration Set of 5 Samples (20g each) of an Allergen in a Specific Matrix is provided for qualitative determination as well as a „Blank Sample“ (Matrix). The allergen levels of the calibration row range from 1/10 to at least twice the Action Level 1 (VITAL concept 3.0) and judgment value of the German Commission ALTS/ALS. Valuation of PT-results is qualitative in Scores of 1-5 (Score 3 = Action Level successfully identified). Quantitative results will be reported for information. The analytical methods are optional.</i>				
<input type="checkbox"/>	ptALM1 (Lot 23Cr)	ALM-Verification Crustaceae : 5 Samples containing King Prawns in Vegan Spread (levels: 10 / 100 / 250 / 375 / 500 mg/kg)	October	400 €
<input type="checkbox"/>	ptALM2 (Lot 23Fi)	ALM-Verification Fish : 5 Samples containing Cod Fish in Wasabi Peas Snack (powder) (levels: 1,3 / 6,5 / 13 / 35 / 70mg/kg)	May	400 €
<i>More Action Level Single-PT-Sets are available for method validation and verification! See also our SPT programme</i>				
<i>Other Elicitors of Allergies and Intolerance Reactions</i>				
<input type="checkbox"/>	ptAI01 (Lot 23LaFr)	Lactose + Fructose in "lactose-free" Food (Baking Mix for Chocolate Cake) (2 different food matrix samples of 25g each, as well as a Spiking Level Sample (25g) with similar levels in simple matrix).	May	280 €

Subscribe	Art.-No. (Lot-No.)	Proficiency Test (PT) 2023	Shipping Samples	Price
Allergen Swab Test (qualitative) - Wipe Test for Surfaces <i>ptASW1-3: There are 2 samples of a Test Surface with Allergens offered for the detection of each allergen by means of a wipe test (swab test). The allergens to be tested were attached to one test surface, while the other test surface is "negative". The evaluation is exclusively qualitative (positive/negative) with indication of the obtained agreement with the consensus values of the participants and the spiking of the surfaces. As an analytical test method, the Lateral-Flow-Test is particularly suitable. ELISA, PCR and LC/MS can also be used.</i>				
<input type="checkbox"/>	ptASW1 (Lot 23)	Allergen Swab Test I: Gluten, Peanut, Milk (incl. β-lactoglobulin, Casein) and Soya (2 Test Surfaces per Allergen, qualitative)	April	360 €
<input type="checkbox"/>	ptASW2 (Lot 23)	Allergen Swab Test II: Crustaceae, Egg, Fish and Mustard (2 Test Surfaces per Allergen, qualitative)	June	360 €
<input type="checkbox"/> NEW	ptASW3 (Lot 23)	Allergen Swab Test III: Almond, Hazelnut, Molluscs and Sesame (2 Test Surfaces per Allergen, qualitative)	July	360 €
Allergen Rinsing Water Test (qualitative + quantitative) <i>ptARW: 3 washing solutions with different allergen contents are offered.</i>				
<input type="checkbox"/>	ptARW1 (Lot 23)	Allergen Rinsing Water Test I: Egg, Gluten, Milk, Mustard and Peanut (3 test solutions with decreasing concentrations, qualitative + quantitative)	March	360 €

2. Mykotoxins - Accredited Programs -

Subscribe	Art.-No. (Lot-No.)	Proficiency Test (PT) 2023	Shipping Samples	Price
All Mycotoxin-PTs (qualitative + quantitative) <ul style="list-style-type: none"> ▶ 2 different samples with "positive" or "negative" contents ▶ qualitative evaluation: agreement pos/neg. in relation to EU Maximum Levels (EU/VO 1881/2006) ▶ quantitative evaluation by z-scores 				
<input type="checkbox"/> up to 3) ¹ <input type="checkbox"/> up to 5) ²	ptMYS1 (Lot 23Multi)	Mycotoxin-Screening: Aflatoxins, Ochratoxin A, Deoxynivalenol, Zearalenone and Fumonisin in Breakfast Cereals (Muesli with corn and/or other cereal products, almonds and dried fruits), 2 different samples each 200g	April	380 € ¹ 480 € ²
<i>1) Determination of up to 3 Parameters, 2) Determination of up to 5 Parameters</i>				
<input type="checkbox"/>	ptMY10 (Lot 23Afl)	Aflatoxin in Nut Butter 2 different samples each 50g	June	290 €
<input type="checkbox"/>	ptMY80 (Lot 23Erg)	Ergot-Alkaloides in Cereals 2 different samples each 80g	December	290 €

3. Potentially toxic Compounds, Contaminants and Residues

Subscribe	Art.-No. (Lot-No.)	Proficiency Test (PT) 2023	Shipping Samples	Price
<input type="checkbox"/>	ptTX01 (Lot 23Cum)	Coumarin in Cinnamon Powder 2 samples each 50g	November	275 €
<input type="checkbox"/>	ptRE01 (Lot 23ClEth)	Ethylene oxide/2-Chloroethanol in Sesame 2 samples each 30g	April	290 €
<input type="checkbox"/>	ptRE02 (Lot 23ClEth)	Ethylene oxide/2-Chloroethanol in Spice Mixture 2 samples each 30g	September	290 €
<input type="checkbox"/>	ptRE01 + ptRE02	Order of both Ethylene Oxide/2-Chloroethanol Proficiency Tests (with special discount)	April + September	480 €

4. Genetically modified Organisms

Subscribe	Art.-No. (Lot-No.)	Proficiency Test (PT) 2023	Shipping Samples	Price
<input type="checkbox"/>	ptGMS1 (Lot 23MaSo)	GMO-Screening I (qualitative): 5 Samples with positive/negative amounts of p-35S, t-NOS, p-FMV, CP4 EPSPS, PAT, Cry1Ab/Ac / GMO-Maize (Bt11, MIR604) and GMO-Soya (RR GTS 40-3-2, RR2 MON89788) (5 different samples 10g each)	May	330 €
<input type="checkbox"/>	ptGMS2 (Lot 23PoRaSu)	GMO-Screening II (qualitative): 5 Samples with positive/negative amounts of p-FMV, CP4 EPSPS, nptII, PAT, BAR / GMO-Potato Amadea (AM04-1020), GMO-Rape Seed (MON88302) / Canola and GMO-Sugar Beet (5 different samples 10g each)	October	330 €
<input type="checkbox"/>	ptGMF (Lot 23MaSoRa)	GMO-Determination in Feed (qualitative + quantitative): GMO-Soya (RR and RR2), GMO-Maize (2 Events), GMO-Rape Seed / Canola (1 Event) 2 different samples 10g each	September	390 €

5. Nutrients, Additives and Authenticity

Subscribe	Art.-No. (Lot-No.)	Proficiency Test (PT) 2023	Shipping Samples	Price
<input type="checkbox"/>	ptAU01 (Lot 23Filn)	Dietary Fiber, soluble / insoluble Fiber and Inulin in Cereal Product (2 different samples, 40g each with and without Inulin)	May	315 €
<input type="checkbox"/>	ptAU05 (Lot 23Stev)	Steviosides in Beverage Powder (2 samples 25g each)	June	250 €
<input type="checkbox"/>	ptAU09 (Lot 23Nic)	Nicotine in E-Cigarette-Liquid (2 samples 10g each)	February	220 €
<input type="checkbox"/> NEW	ptAU11 (Lot 23Veg)	Vegan Compliance Check – 2 Samples qualitative: Animal ingredients in Vegan Burger (fried) Negative sample and positive sample with at least two animal components > 0.1% (e.g. egg, fish, meat/gelatine, milk components and/or honey) (2 different samples of 50 g each)	April	350 €
Animal Species-Screening (qualitative) <i>ptAUS1-5: 4 Samples or 6 Samples (ptAUS4) with portions of Animal Origin (meat, milk, fish, crustacean or mollusc products) are provided. The evaluation of results is strictly qualitative with indication of the obtained agreements with the consensus values of the participants and the spiking of samples. The analytical methods are optional.</i>				
<input type="checkbox"/>	ptAUS1 (Lot 23Meat)	Animal Species-Screening I – 4 Samples qualitative: Beef, Poultry (Duck (genus: <i>anas</i> or <i>cairina</i>) and Chicken), Sheep and Goat in cooked Meat Product (Sausage Meat, 4 different samples 30g each)	March	360 €
<input type="checkbox"/>	ptAUS2 (Lot 23Meat)	Animal Species-Screening II – 4 Samples qualitative: Donkey, Horse, Game (Deer (<i>Capreolus capreolus</i>) or fallow deer (<i>Dama dama</i>), Pork and Poultry (Turkey) in cooked Meat Product (Sausage Meat, 4 different samples 30g each)	October	360 €
<input type="checkbox"/>	ptAUS3 (Lot 23Milk)	Animal Species-Screening III – 4 Samples qualitative: Buffalo, Cow's, Sheep's and Goat's Milk in Dairy Product (Mozzarella Cheese) (freeze-dried Mixtures, 4 different Samples 20g each)	June	380 €
<input type="checkbox"/>	ptAUS4 (Lot 23SF)	Seafood-Screening – 6 Samples qualitative: Bonito (<i>Katsuwonus</i> <i>pelamis</i>), Cod Fish (<i>Gadus morhua</i>), Rainbow Trout (<i>Oncorhynchus</i> <i>mykiss</i>), White Tiger Shrimp (<i>Penaeus vannamei</i>) and one of the Cephalopods: Octopus (<i>Octopus vulgaris</i>) or Jumbo Squid (<i>Dosidicus</i> <i>gigas</i>) freeze-dried, 5 samples (10g each) with only one of the indicated species, 6th sample (10g) mixed sample of 2-4 species	November	410 €
<input type="checkbox"/>	ptAUS5 (Lot 23PF)	Petfood-Screening III – 4 Samples qualitative: Elk, Lamb, Ostrich, Rabbit and Beef in Dry Dog Food (4 different samples 30g each)	February	360 €

6. Food Supplements and Dietetic Products

Subscribe	Art.-No. (Lot-No.)	Proficiency Test (PT) 2023	Shipping Samples	Price
<input type="checkbox"/>	ptSU01 (Lot 23WaVi)	Food Supplement I: Vitamins B1, B2, B6, B12, Biotin, Vitamin C, Folic Acid, Niacin and Pantothenic Acid (enriched / high levels, 2 samples 50g each)	March	370 €
<input type="checkbox"/>	ptSU03 (Lot 23Min)	Food Supplement II: B, Ca, Cr, Cu, Fe, K, Mg, Mn, Mo, P, Se, Zn (enriched / low levels, 2 samples 50g each)	May	370 €
<input type="checkbox"/>	ptSU07 (Lot 23FaVi)	Dietetic Food I: Vitamine A, E, D3, K1 and β-Carotene in Meal Replacement (enriched / low levels, 2 samples 50g each)	October	370 €
<input type="checkbox"/>	ptSU08 (Lot 23Elem)	Minerals, Trace Elements and Heavy Metals in Infant Food , approx. 20 Elements incl. As, Cd, Hg and Pb (mostly "natural" / partly low levels, 2 samples 8g each)	September	370 €
<input type="checkbox"/>	ptSU10 (Lot 23FrAa)	Free Amino Acids and Taurine in Food (Infant Food) (2 samples 10g each)	September	390 €

7. Cosmetic Products and Contact Materials

Subscribe	Art.-No. (Lot-No.)	Proficiency Test (PT) 2023	Shipping Samples	Price
<input type="checkbox"/>	ptCO01 (Lot 23ColPre)	Cosmetic Products I: Coffein, Panthenol and Preservatives (Benzoic Acid, Methyl-, Propylparaben, Phenoxyethanol, Sorbic Acid) in Hair Shampoo (2 samples 50g each)	April	320 €
<input type="checkbox"/>	ptCO03 (Lot 23Vit)	Cosmetic Products III: Coenzym Q10, Panthenol, Retinol and Tocopherol in Skin Cream (2 samples 25g each)	November	280 €
<input type="checkbox"/>	ptCO06 (Lot 23UV)	Cosmetic Products II: UV-Filter in Sun Protection Product Spray (Liquid solution, 2 samples 25mL each)	February	280 €
<input type="checkbox"/>	ptCM05 (Lot 23Elem)	Children's Product I: Analysis of 8 Elements in Make-Up (As, Ba, Cd, Cr, Hg, Pb, Sb, Se) 2 samples 10g each	May	320 €
<input type="checkbox"/>	<i>Please indicate!</i>	Double Sample Material available! <i>You can now order twice the sample quantity for all our PTs. Please enter the requested DLA-No. below:</i>	<i>Mit EP-Probenversand</i>	90 €
<i>After the evaluation report is finished our PT materials are usually available as Reference Material! More information on our website!</i>				
<i>More parameter/matrices are available in Single-Proficiency Tests! See also our SPT programme on our website</i>				

All prizes plus tax and plus shipping costs.

Our general terms and conditions apply, see at www.dla-lvu.de

Shipment currently: € 12,50 in Germany, € 20,00 into EU-Countries and Europe, € 28,00 North America
 (Express delivery on request plus € 35,00), € 80,00 other countries (for example in South America and Asia)

Special Offer in 2023:
5% discount for booking ≥ 10 registrations
(valid for PT annual registrations until 15th January 2023)

Contact Information (Please send the marked form including the last page!)

In case the customer number is known the address is not required.

After receiving your registration, we will send you a registration confirmation by e-mail.

Institution / Company	<input type="text"/>
Customer number * (If known)	<input type="text"/>
VAT number	<input type="text"/>
Street / Post Box	<input type="text"/>
Code City Country	<input type="text"/>
Order Information	<input type="text"/>
	(Please fill in your internal order information and/or how did you hear about us!)
Invoice Address (if different from above)	<input type="text"/>
Invoices via E-Mail to	<input type="text"/>
Telephone / Telefax	<input type="text"/>
E-Mail	<input type="text"/>
Contact Person	<input type="text"/>
*	(in case DLA knows your further contact details only the Customer Number is required!)
<input type="checkbox"/>	I consent to the Data Protection Policy on www.dla-lvu.de
<input type="checkbox"/>	Hereby we subscribe to the marked proficiency tests of the DLA programs 2023. Cancellation is possible free of charge up to 2 months before shipment.

Registration Form: DLA - Proficiency Tests 2023

Please send the PDF-form by email to pt@dla-lvu.de
or register directly using our website form